

SVEMA SANITRUST WIPES






Surface Sanitising Wipes White & Surface Sanitising Wipes Blue

Ethanol-based rapid disinfection wipes for alcohol-tolerant surfaces in food, veterinary, medical, industrial, domestic and institutional areas.

Sanitrust Wipes have a broad spectrum of bactericidal (including MRSA, VRE), mycobactericidal, fully virucidal (incl. all enveloped viruses, Noro, Rota, Adeno, Polio), fungicidal and yeasticidal activity.



TECHNICAL INFORMATION

Active substance	<ul style="list-style-type: none"> Ethanol CAS no 64-17-5 (57 g); Propan-2-ol CAS no 67-63-0 (6 g) 				
Specifications	<ul style="list-style-type: none"> Fragrance: unfragranced Canister of 150 pcs blue wipes: 20x20cm; 30g/m²; nonwoven, lint-free; Canister of 200 pcs white wipes: 20x20cm; 28g/m²; nonwoven, lint-free 				
Preservation	<ul style="list-style-type: none"> Can be kept for 30 months in its original closed packaging under normal storage conditions. Do not freeze, store between 5°C and 30°C 				
Directions for use	1. To pull the first wipe, remove the canister lid and pull the first wipe from the centre of the reel with clean hands. Replace the canister lid with the first wipe pulled through the opening in the lid. 	2. Wipe pre-cleaned surface. Leave it to act 5 minutes. Close the lid after use. 	3. Leave to dry, do not rinse. Use several wipes if necessary. 	4. Dispose of wipes according to the instructions on the product label. 	5. Repeat the operation in accordance with the hygiene protocol in place. 
Regulation	<ul style="list-style-type: none"> Biocide PT2, PT4. AL - other liquid to be applied undiluted Professional use 				
Compatibility	<ul style="list-style-type: none"> The product is compatible with most smooth surfaces except alcohol-sensitive surfaces (plexiglass, etc.) The product, in its use as a food contact surface sanitiser, does not alter the taste or smell of food, according to Tetrad Test ASTM International E3009-15 				
Skin compatibility	<ul style="list-style-type: none"> Good skin compatibility Tested under dermatological control 				
Precaution for use	<ul style="list-style-type: none"> Do not swallow Do not mix with other products Use biocides safely Always read the label and product Information before use Dispose of contents/container into an approved waste disposal container Do not reuse packaging For further information, refer to the product safety data sheet In the event of an emergency, contact the nearest Poison Centre 				

ANTIMICROBIAL ACTIVITY	
Active against	Contact time
BACTERIA (INCL. MRSA, VRE):	15 sec
YEASTICIDAL ACTIVITY: CANDIDA ALBICANS	15 sec
MYCOBACTERIA: (M.AVIUM, M.TERRAE)	30 sec
FUNGICIDAL ACTIVITY: ASBERGILLUS BRASILIENSIS	120 sec
VIRUSES: VACCINIA, HIV, HBV, HCV, ROTA, NORO	15 sec
VIRUSES: ADENOVIRUS	30 sec
VIRUSES: POLIO VIRUS	120 sec
STANDARDS	EN13727, EN14561, EN13697, EN13624, EN14562, EN14348, EN14563, EN14476, EN16777, EN17387, EN16615, EN1650, EN1276.

*Study conducted in clean conditions

PACKAGING	SKU CODE	CASE
Product name		Pieces
SVEMA SANITRUST WIPES (200 Count Canister)	5020-80-06	6
SVEMA SANITRUST WIPES BLUE (150 Count Canister)	5030-80-06	6